

ENDOTRACHEAL TUBE

ANAESTHESIA

- Safe & effective airway management.
- Suitable for general & emergency applications.
- User friendly for anaesthetist & patients.
- Easy insertion and less trauma.

- **Smoothly connected inflation tube** maximizes patient's comfort.

- Smooth, flexible and kink resistant inflation line ensures comfort and **effective inflation & deflation process.**
- Size marking on pilot balloon for better **product identification.**



- An additional opening to the tip, the Murphy eye has a polished smoothness ensuring atraumatic intubation and extubation.


- Radio-opaque line allowing **clear identification** of the tube placement on radiographic images.

- Thermosensitive tubing softens at body temperature to conform to the anatomy of the respiratory tract.

- Double-ring marks help guide intubation through laryngoscopy.

- Smooth bull-nose tip enables **atraumatic intubation.** Short tip minimizes contact with tracheal mucosa and increases the distance from the carina.

- Soft, thin-walled cuff for effective sealing. Heat-welded cuff ensures **smooth cuff-to-tube transition.** For atraumatic intubation and extubation.



WIRE REINFORCED ENDOTRACHEAL
TUBE CUFFED



WIRE REINFORCED ENDOTRACHEAL TUBE
CUFFED WITH STYLET

The volume, low pressure cuff, soft, low trauma tip



Printed pilot balloon for better traceability



The integrated stainless steel reinforcing wire provides a safe, kink resistant airway regardless of head and neck position.



Ideal Healthcare Sdn. Bhd.(484861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 422 7775 / 422 8775 Fax : 6 (04) - 422 6775
Email : info@idealcare.com.my

E. L. NO.: KP51039514




www.idealcare.com.my

0120

idealcareMY

IHSB-CAT-BDD-TDR-ANS-05

 Regular HVLP Cuffed	Catalogue No.	Size	Catalogue No.	Size	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61110030	3.0 mm	IC61110070	7.0 mm	
	IC61110035	3.5 mm	IC61110075	7.5 mm	
	IC61110040	4.0 mm	IC61110080	8.0 mm	
	IC61110045	4.5 mm	IC61110085	8.5 mm	
	IC61110050	5.0 mm	IC61110090	9.0 mm	
	IC61110055	5.5 mm	IC61110095	9.5 mm	
	IC61110060	6.0 mm	IC61110100	10.0 mm	
 Regular SVLP Cuffed	IC61110065	6.5 mm	-	-	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61150045	4.5 mm	IC61150070	7.0 mm	
	IC61150050	5.0 mm	IC61150075	7.5 mm	
	IC61150055	5.5 mm	IC61150080	8.0 mm	
	IC61150060	6.0 mm	IC61150085	8.5 mm	
 Wired Reinforced HVLP Cuffed	IC61150065	6.5 mm	IC61150090	9.0 mm	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61214030	3.0 mm	IC61214065	6.5 mm	
	IC61214035	3.5 mm	IC61214070	7.0 mm	
	IC61214040	4.0 mm	IC61214075	7.5 mm	
	IC61214045	4.5 mm	IC61214080	8.0 mm	
	IC61214050	5.0 mm	IC61214085	8.5 mm	
	IC61214055	5.5 mm	IC61214090	9.0 mm	
 Wired Reinforced HVLP Cuffed With Stylet	IC61214060	6.0 mm	IC61214095	9.5 mm	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61214130	3.0 mm	IC61214165	6.5 mm	
	IC61214135	3.5 mm	IC61214170	7.0 mm	
	IC61214140	4.0 mm	IC61214175	7.5 mm	
	IC61214145	4.5 mm	IC61214180	8.0 mm	
	IC61214150	5.0 mm	IC61214185	8.5 mm	
	IC61214155	5.5 mm	IC61214190	9.0 mm	
 Oral Curved SVLP Cuffed	IC61214160	6.0 mm	IC61214195	9.5 mm	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61213040	4.0 mm	IC61213070	7.0 mm	
	IC61213045	4.5 mm	IC61213075	7.5 mm	
	IC61213050	5.0 mm	IC61213080	8.0 mm	
	IC61213055	5.5 mm	IC61213085	8.5 mm	
	IC61213060	6.0 mm	IC61213090	9.0 mm	
 Nasal Curved SVLP Cuffed	IC61213065	6.5 mm	-	-	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61313040	4.0 mm	IC61313070	7.0 mm	
	IC61313045	4.5 mm	IC61313075	7.5 mm	
	IC61313050	5.0 mm	IC61313080	8.0 mm	
	IC61313055	5.5 mm	IC61313085	8.5 mm	
	IC61313060	6.0 mm	IC61313090	9.0 mm	
 Regular Plain Uncuff	IC61313065	6.5 mm	-	-	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61130020	2.0 mm	IC61130065	6.5 mm	
	IC61130025	2.5 mm	IC61130070	7.0 mm	
	IC61130030	3.0 mm	IC61130075	7.5 mm	
	IC61130035	3.5 mm	IC61130080	8.0 mm	
	IC61130040	4.0 mm	IC61130085	8.5 mm	
	IC61130045	4.5 mm	IC61130090	9.0 mm	
	IC61130050	5.0 mm	IC61133095	9.5 mm	
	IC61130055	5.5 mm	IC61130100	10.0 mm	
 Oral Curved Plain Uncuff	IC61130060	6.0 mm	-	-	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61233030	3.0 mm	IC61233065	6.5 mm	
	IC61233035	3.5 mm	IC61233070	7.0 mm	
	IC61233040	4.0 mm	IC61233075	7.5 mm	
	IC61233045	4.5 mm	IC61233080	8.0 mm	
	IC61233050	5.0 mm	IC61233085	8.5 mm	
	IC61233055	5.5 mm	IC61233090	9.0 mm	
 Nasal Curved Plain Uncuff	IC61233060	6.0 mm	IC61233090	9.0 mm	Packaging: 10 pcs / box, sterile Shipping carton: 10 boxes / carton
	IC61333030	3.0 mm	IC61333065	6.5 mm	
	IC61333035	3.5 mm	IC61333070	7.0 mm	
	IC61333040	4.0 mm	IC61333075	7.5 mm	
	IC61333045	4.5 mm	IC61333080	8.0 mm	
	IC61333050	5.0 mm	IC61333085	8.5 mm	
 Intubation Stylet	IC61333055	5.5 mm	IC61333090	9.0 mm	Packaging: 20 pcs / box, sterile, 20 boxes / carton
	IC61333060	6.0 mm	-	-	
	IC64200006	6Fr (2.0 mm)	Length: 305 mm		
	IC64200010	10Fr (3.3 mm)	Length: 390 mm		
	IC64200014	14Fr (4.7 mm)	Length: 420 mm		