


MEDICAL FACE SHIELD / ISOLATION MASK

PPE COVID 19

- Provide barrier protection to general public and healthcare personnel as PPE against splashes and spatter of body fluids from contacting directly with the wearer.
- Effectively prevents from saliva splashing and fully protects the face and eyes.
- Anti-fog which effectively prevent fog and maintain good visibility.
- Clear face shield film.
- Elastic band and comfortable foam.
- Super light and thin 18mm thickness.
- Size: 33cm x 22cm.
- Single use.
- Non-sterile.



- 1 Clean face shield film.
- 2 Foam.
- 3 Elastic band.
- 4 Bonding point.
- 5 Shield label.

IHSB REF. NO.	SIZE
IC0032	22cm x 33cm



Ideal Healthcare Sdn. Bhd.(484861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 422 7775 / 422 8775 Fax : 6 (04) - 422 6775
Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

"Your One Stop Centre for Medical & Healthcare Products"

www.idealcare.com.my



idealcareMY

IHSB-CAT-BDD-PPE-10-00