

IDEALCARE®
.....the quality you can trust

PRODUCT CATALOGUE
2022/2023



Scan Me



catalogue

Catalogue Version 4.0

CONTENTS

MANUFACTURER PRODUCTS

IDEALCARE®
.....the quality you can trust

• Anaesthesia	1-4
• Respiratory Care	4-5
• Urology	6-8
• Infusion & Transfusion	8-12
• Gastroenterology	12
• Surgery & Wound Drainage	13-14
• Pathology	14-15
• Obstetrics & Gynaecology	16
• Procedure & Clinical Sets	16-18
• Procedure Set Components	18
• Wound Care & Dressing	19-22
• Pharmaceutical Plastic Container	22-23
• Rehabilitation	23
• Others	23
• PPE Covid19	24-27

AUTHORIZE REPRESENTATIVE

POLYMED
Medical Devices

• Anaesthesia	28
• Infusion & Transfusion	27-31
• Others	31
• Dialysis	32
• Gastroenterology	32
• Surgery & Wound Drainage	32
• Central Venous Access Catheter	33

Endotracheal Tube (MDA. Reg. No: GB49099443317)


 Endotracheal Tube Wire
Reinforced Uncuff

 Wire Reinforced Endotracheal Tubes
Cuffed

 Wire Reinforced Endotracheal Tubes
Cuffed With Stylet

	CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE	SKU
 REGULAR HVLP CUFFED	IC 61110030	3.0 mm	IC 61110065	6.5 mm	Packaging: 10 pcs / box, sterile
	IC 61110035	3.5 mm	IC 61110070	7.0 mm	
	IC 61110040	4.0 mm	IC 61110075	7.5 mm	
	IC 61110045	4.5 mm	IC 61110080	8.0 mm	
	IC 61110050	5.0 mm	IC 61110085	8.5 mm	
	IC 61110055	5.5 mm	IC 61110090	9.0 mm	
	IC 61110060	6.0 mm			
 REGULAR SVLP CUFFED	IC 61150045	4.5 mm	IC 61150070	7.0 mm	Packaging: 10 pcs / box, sterile
	IC 61150050	5.0 mm	IC 61150075	7.5 mm	
	IC 61150055	5.5 mm	IC 61150080	8.0 mm	
	IC 61150060	6.0 mm	IC 61150085	8.5 mm	
	IC 61150065	6.5 mm	IC 61150090	9.0 mm	
 WIRED REINFORCED HVLP CUFFED	IC 61214030	3.0 mm	IC 61214065	6.5 mm	Packaging: Single pack, sterile
	IC 61214035	3.5 mm	IC 61214070	7.0 mm	
	IC 61214040	4.0 mm	IC 61214075	7.5 mm	
	IC 61214045	4.5 mm	IC 61214080	8.0 mm	
	IC 61214050	5.0 mm	IC 61214085	8.5 mm	
	IC 61214055	5.5 mm	IC 61214090	9.0 mm	
	IC 61214060	6.0 mm	IC 61214095	9.5 mm	
 WIRED REINFORCED HVLP CUFFED WITH STYLET	IC61214130	3.0 mm	IC61214165	6.5 mm	Packaging: Single pack, sterile
	IC61214135	3.5 mm	IC61214170	7.0 mm	
	IC61214140	4.0 mm	IC61214175	7.5 mm	
	IC61214145	4.5 mm	IC61214180	8.0 mm	
	IC61214150	5.0 mm	IC61214185	8.5 mm	
	IC61214155	5.5 mm	IC61214190	9.0 mm	
	IC61214160	6.0 mm	IC61214195	9.5 mm	
 ORAL CURVED SVLP CUFFED	IC61213040	4.0 mm	IC61213070	7.0 mm	Packaging: Single pack, sterile
	IC61213045	4.5 mm	IC61213075	7.5 mm	
	IC61213050	5.0 mm	IC61213080	8.0 mm	
	IC61213055	5.5 mm	IC61213085	8.5 mm	
	IC61213060	6.0 mm	IC61213090	9.0 mm	
	IC61213065	6.5 mm			

	CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE	SKU
 NASAL CURVED SVLP CUFFED	IC61313045	4.5 mm	IC61313070	7.0 mm	Packaging: Single pack, sterile
	IC61313050	5.0 mm	IC61313075	7.5 mm	
	IC61313055	5.5 mm	IC61313080	8.0 mm	
	IC61313060	6.0 mm	IC61313085	8.5 mm	
	IC61313065	6.5 mm	IC61313090	9.0 mm	
 WIRED REINFORCED PLAIN UNCUFF	IC 61134030	3.0 mm	IC 61134065	6.5 mm	Packaging: Single pack, sterile
	IC 61134035	3.5 mm	IC 61134070	7.0 mm	
	IC 61134040	4.0 mm	IC 61134075	7.5 mm	
	IC 61134045	4.5 mm	IC 61134080	8.0 mm	
	IC 61134050	5.0 mm	IC 61134085	8.5 mm	
	IC 61134055	5.5 mm	IC 61134090	9.0 mm	
 REGULAR PLAIN UNCUFF	IC 61130020	2.0 mm	IC 61130065	6.5 mm	Packaging: 10 pcs / box, sterile
	IC 61130025	2.5 mm	IC 61130070	7.0 mm	
	IC 61130030	3.0 mm	IC 61130075	7.5 mm	
	IC 61130035	3.5 mm	IC 61130180	8.0 mm	
	IC 61130040	4.0 mm	IC 61130185	8.5 mm	
	IC 61130045	4.5 mm	IC 61130190	9.0 mm	
	IC 61130050	5.0 mm	IC 61130195	9.5 mm	
	IC 61130055	5.5 mm	IC 61130100	10.0 mm	
	IC 61130060	6.0 mm			
	 ORAL CURVED PLAIN UNCUFF	IC61233030	3.0 mm	IC61233065	
IC61233035		3.5 mm	IC61233070	7.0 mm	
IC61233040		4.0 mm	IC61233075	7.5 mm	
IC61233045		4.5 mm	IC61233080	8.0 mm	
IC61233050		5.0 mm	IC61233085	8.5 mm	
IC61233055		5.5 mm	IC61233090	9.0 mm	
IC61233060		6.0 mm			
 NASAL CURVED PLAIN UNCUFF	IC61333030	3.0 mm	IC61333065	6.5 mm	Packaging: Single pack, sterile
	IC61333035	3.5 mm	IC61333070	7.0 mm	
	IC61333040	4.0 mm	IC61333075	7.5 mm	
	IC61333045	4.5 mm	IC61333080	8.0 mm	
	IC61333050	5.0 mm	IC61333085	8.5 mm	
	IC61333055	5.5 mm	IC61333090	9.0 mm	
	IC61333060	6.0 mm			
 INTUBATION STYLET	IC64200006	6Fr (2.0 mm)	Packaging: Single pack, sterile		
	IC64200010	10Fr (3.3 mm)			
	IC64200014	14Fr (4.7 mm)			

Tracheostomy Tube (MDA. Reg. No: GB8479366216)



Cuffed

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC62520045	4.5 mm	IC62520075	7.5 mm
IC62520050	5.0 mm	IC62520080	8.0 mm
IC62520055	5.5 mm	IC62520085	8.5 mm
IC62520060	6.0 mm	IC62520090	9.0 mm
IC62520065	6.5 mm	IC62520095	9.5 mm
IC62520070	7.0 mm	IC62520100	10.0 mm
Packaging: Single pack, sterile			



Plain Uncuff

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC62540030	3.0 mm	IC62540070	7.0 mm
IC62540035	3.5 mm	IC62540075	7.5 mm
IC62540040	4.0 mm	IC62540080	8.0 mm
IC62540045	4.5 mm	IC62540085	8.5 mm
IC62540050	5.0 mm	IC62540090	9.0 mm
IC62540055	5.5 mm	IC62540095	9.5 mm
IC62540060	6.0 mm	IC62540100	10.0 mm
IC62540065	6.5 mm	Packaging: Single pack, sterile	

Laryngeal Airway Mask

(MDA. Reg. No: GB67609671618)



Disposable (PVC)		
CATALOGUE NO.	SIZE	PATIENT WEIGHT
IC6310D	1.0 mm	Neonate < 5 kgs
IC6315D	1.5 mm	Infant, 5-10 kgs
IC6320D	2.0 mm	Child, 10-20 kgs
IC6325D	2.5 mm	Child, 20-30 kgs
IC6330D	3.0 mm	Child, 30-50 kgs
IC6340D	4.0 mm	Adult, 50-70 kgs
IC6350D	5.0 mm	Adult, 70-100 kgs
Packaging: Single pack, sterile		



Reusable (Silicone)		
CATALOGUE NO.	SIZE	PATIENT WEIGHT
IC6410R	1.0mm	Neonate < 5 kgs
IC6415R	1.5mm	Infant, 5-10 kgs
IC6420R	2.0mm	Child, 10-20 kgs
IC6425R	2.5mm	Child, 20-30 Kgs
IC6430R	3.0mm	Child, 30-50 Kgs
IC6440R	4.0mm	Adult, 50-70 kgs
IC6450R	5.0mm	Adult, 70-100 kgs
Packaging: Single pack, sterile		

Suction Catheter with Finger Tip Control

(MDA. Reg. No: GB14585265317)



CATALOGUE NO.	SIZE	COLOR CODE
IC2030-06	6Fr	Green
IC2030-08	8Fr	Blue
IC2030-10	10Fr	Black
IC2030-12	12Fr	White
IC2030-14	14Fr	Green
IC2030-16	16Fr	Orange
IC2030-18	18Fr	Red

Packaging: 100 pcs / box, sterile

RESPIRATORY CARE

Silicone Resuscitator

(MDA. Reg. No: GB181681028318)



CATALOGUE NO.	SIZE
IC2103	Adult
IC2203	Child
IC2303	Infant

Packaging: Single pack

Nasal Oxygen Cannula

(MDA. Reg. No: GB37077264717)



CATALOGUE NO.	SIZE
IC2016A	Adult
IC2016P	Paediatric
IC2016I	Infant
IC2016N	Neonate

Packaging: Single pack

(MDA. Reg. No: GB422641130318)

CATALOGUE NO.	SIZE
IC2017A	Adult
IC2017P	Paediatric

Packaging: Single pack, sterile

Oxygen Reservoir

(MDA. Reg. No: GMD39866449218A)



CATALOGUE NO.	SIZE
IC0438	2500ml
IC0440	6000ml

Packaging: Single pack

Oxygen Mask with Tubing (MDA. Reg. No: GB17840442917)


CATALOGUE NO.	SIZE
IC2020A	Elongated Adult
IC2020P	Elongated Paediatric
IC2020 I	Standard Infant
IC2020N	Standard Neonate

Packaging: Single pack

Oxygen Mask High Concentration Non Rebreathing Mask

(MDA. Reg. No: GB828721028518)



CATALOGUE NO.	SIZE
IC1110A	Adult
IC1112P	Paediatric

Packaging: Single pack

Nebulizer Mask with Chamber & Tubing (MDA. Reg. No: GB96047264817)


CATALOGUE NO.	SIZE
IC1120	Elongated Adult
IC1122	Elongated Paediatric
IC1123	Standard Neonate

Packaging: Single pack

Oxygen Air Entrainment (MDA. Reg. No: GB195711256919)


CATALOGUE NO.	SIZE
IC7411	Adult
IC7412	Paediatric

Packaging: Single pack

Adjusted Oxygen Concentration	Color Code
24% (Blue)	
26% (Yellow)	
35% (Dark Green)	
40% (Pink)	
60% (Red)	

Heat Moisture Exchange (HME) Filter (MDA. Reg. No: GMD4748269616A)


CATALOGUE NO.	SIZE
IC158000	Standard

Packaging: Single pack, sterile

Latex Foley Catheter, Sterile

(MDA. Reg. No: GB9868069116)

All Silicone Foley Catheter, Sterile

(MDA. Reg. No: GC8986604218)



2-Way Standard Latex Foley Catheter Silicone Coated

2-Way Premium Gold Foley Catheter

2-Way All Silicone Foley Catheter

2-Way Paediatric Latex Foley Catheter Silicone Coated	2-Way Standard Latex Foley Catheter Silicone Coated			3-Way Standard Latex Foley Catheter Silicone Coated			Packaging: 10 pcs / box, sterile
2-Way Paediatric (HV)	Note: HV = Hard Valve						
	2-Way Standard (HV)			3-Way Standard (HV)			
3-5ml/cc	5-15ml/cc	20-30ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc	Size (Ch/Fr)	
IC23120605	-	-	-	-	-	06	
IC23120805	-	-	-	-	-	08	
IC23121005	-	-	-	-	-	10	
-	IC21111210	IC21111230	-	-	-	12	
-	IC21111410	IC21111430	-	-	-	14	
-	IC21111610	-	IC21111830	IC25121610	IC25121630	16	
-	IC21111810	-	IC21111830	IC25121810	IC25121830	18	
-	IC21112010	-	IC21112030	IC25122010	IC25122030	20	
-	IC21112210	-	IC21112230	IC25122210	IC25122230	22	
-	IC21112410	-	IC21112430	IC25122410	IC25122430	24	
-	IC21112610	-	IC21112630	IC25122610	IC25122630	26	
-	IC21112810	-	IC21112830	-	-	28	
-	-	-	IC21113030	-	-	-	

2-Way Premium Paediatric Gold Foley Catheter	2-Way Premium Gold Foley Catheter		3-Way Premium Gold Foley Catheter			Packaging: 10 pcs / box, sterile
3-5ml/cc	5-15ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc	Size (Ch/Fr)	
GD23120605	-	-	-	-	06	
GD23120805	-	-	-	-	08	
GD23121005	-	-	-	-	10	
-	GD21111210	GD21111230	-	-	12	
-	GD21111410	GD21111430	-	-	14	
-	GD21111610	GD21111630	GD25121610	GD25121630	16	
-	GD21111810	GD21111830	GD25121810	GD25121830	18	
-	GD21112010	GD21112030	-	GD25122030	20	
-	GD21112210	GD21112230	GD25122210	GD25122230	22	
-	GD21112410	GD21112430	GD25122410	GD25122430	24	
-	GD21112610	GD21112630	GD25122610	GD25122630	26	

2-Way Paediatric All Silicone Foley Catheter		2-Way All Silicone Foley Catheter		3-Way All Silicone Foley Catheter			Packaging: Single pack, sterile
Note: HV = Hard Valve		Note: HV = Hard Valve		Note: HV = Hard Valve			
2-Way Paediatric (HV)		2-Way (HV)		3-Way (HV)			
1-5ml/cc	3-5ml/cc	5-15ml/cc	30-45ml/cc	5-15ml/cc	30-45ml/cc	Size (Ch/Fr)	
IC43150603	-	-	-	-	-	06	
-	IC43150803	-	-	-	-	08	
-	IC43151003	-	-	-	-	10	
-	-	IC41151205	-	-	-	12	
-	-	IC41151405	-	-	-	14	
-	-	IC41151605	IC41151630	-	IC45151630	16	
-	-	IC41151805	IC41151830	IC45151805	IC45151830	18	
-	-	IC41152005	IC41152030	IC45152005	IC45152030	20	
-	-	IC41152205	IC41152230	-	IC45152230	22	
-	-	IC41152405	IC41152430	-	IC45152430	24	
-	-	IC41152605	IC41152630	-	IC45152630	26	

2-Way Premium Hematuria Reinforced Foley Catheter				Size (Ch/Fr)	Packaging: 10 pcs / box, sterile
Couveleire		Dufour			
30-50 ml/cc	60-80 ml/cc	30-50 ml/cc	60-80 ml/cc		
GD21421830				18	
GD21422030	GD21422060	GD21522030	GD21522060	20	
GD21422230	GD21422260	GD21522230	GD21522260	22	
GD21422430	GD21422460	GD21522430	GD21522460	24	
GD21422630	GD21422660	GD21522630	GD21522660	26	

3-Way Premium Hematuria Reinforced Foley Catheter				Size (Ch/Fr)	Packaging: 10 pcs / box, sterile
Couveleire		Dufour			
30-50 ml/cc	60-80 ml/cc	30-50 ml/cc	60-80 ml/cc		
GD25421830				18	
GD25422030	GD25422060	GD25522030	GD25522060	20	
GD25422230	GD25422260	GD25522230	GD25522260	22	
GD25422430	GD25422460	GD25522430	GD25522460	24	
GD25422630	GD25422660	GD25522630	GD25522660	26	

Nelaton Premium Foley Catheter		Packaging: 10 pcs / box, sterile
CATALOGUE NO.	SIZE	
GD36101200	12	
GD36101400	14	
GD36101600	16	
GD36101800	18	
GD36102000	20	
GD36102200	22	
GD36102400	24	
GD36102600	26	

Male External Catheter (MDA. Reg. No: GMD16184357417A)



CATALOGUE NO.	SIZE	LENGTH
IC1010-S	Small	25mm
IC1010-M	Medium	30mm
IC1010-L	Large	35mm
IC1010-XL	Extra Large	40mm
Packaging: 50 pcs / box, sterile		

Adult Urine Collecting Bag (MDA. Reg. No: GMD91401113817A)



CATALOGUE NO.	SIZE
IC1064(A)	2000 ml
IC1064SH(A)	2000 ml (with urine bag hanger)
Packaging: Single pack, sterile	

Paediatric Urine Collector (MDA. Reg. No: GMD91401113817A)



CATALOGUE NO.	SIZE
IC1062	100 ml
Packaging: 100 pcs / box, sterile	

Urine Collecting Bag with Measured Volume Meter (Urimeter)

(MDA. Reg. No: GMD54729189217A)



CATALOGUE NO.	BAG CAPACITY	CHAMBER CAPACITY
IC1065	2000 ml	250 ml
Packaging: Single pack, sterile		

Urine Collecting Bag with Measured Volume Meter (Urimeter)

(MDA. Reg. No: GMD54729189217A)



CATALOGUE NO.	BAG CAPACITY	CHAMBER CAPACITY
IC1066	2000 ml	500 ml
Packaging: Single pack, sterile		

Urine Bag Hanger (MDA. Reg. No: GMD7717484818A)



CATALOGUE NO.	COLOR
IC1060	Blue
Packaging: Single pack / pack	

INFUSION & TRANSFUSION

I.V. Catheter for Quick Flashback with Adva Needle Technology

(MDA. Reg. No: GMD7717484818A)



CATALOGUE NO.	SIZE	COLOR CODE	ID (mm)	OD (mm)	LENGTH (cm)	FLOW RATE (ml/min)
IC1039-14	14G x 5.1cm	Orange	1.7	2.0	5.1	305
IC1039-16	16G x 5.0cm	Grey	1.3	1.7	5.0	200
IC1039-18	18G x 4.5cm	Green	0.9	1.3	4.5	95
IC1039-20	20G x 3.2cm	Red	0.8	1.1	3.2	65
IC1039-22	22G x 2.5cm	Blue	0.6	0.9	2.5	36
IC1039-24	24G x 1.9cm	Yellow	0.5	0.7	1.9	23
IC1039-26	26G x 1.9cm	Purple	0.45	0.6	1.9	17
Packaging: 50 pcs / box, sterile						

INFUSION & TRANSFUSION

Disposable Intravenous Catheter With Injection Port & Wings

(MDA. Reg. No: GB24437168617)



CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (cm)	FLOW RATE (ml/min)
IC2001-14	14G	Orange	2.0	1.7	5.0	305
IC2001-16	16G	Grey	1.7	1.3	5.0	200
IC2001-18	18G	Green	1.3	0.9	4.5	95
IC2001-20	20G	Red	1.1	0.8	3.3	65
IC2001-22	22G	Blue	0.9	0.6	2.5	36
IC2001-24	24G	Yellow	0.7	0.5	1.9	23
IC2001-26	26G	Purple	0.6	0.45	1.9	17

Packaging: 50 pcs / box, sterile

Disposable Intravenous Catheter With Injection Port & Wings (improved)

(MDA. Reg. No: GB24437168617)



CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (cm)	FLOW RATE (ml/min)
IC 2001-14	14G	Orange	2.0	1.7	5.0	305
IC 2001-16	16G	Grey	1.7	1.3	5.0	200
IC 2001-18	18G	Green	1.3	0.9	4.5	95
IC 2001-20	20G	Red	1.1	0.8	3.3	65

Packaging: 50 pcs / box, sterile

Disposable Intravenous Catheter With With Injection Valve

(MDA. Reg. No: GB24437168617)



CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (cm)	FLOW RATE (ml/min)
IC2002-14	14G	Orange	2.0	1.7	5.0	305
IC2002-16	16G	Grey	1.7	1.3	5.0	200
IC2002-18	18G	Green	1.3	0.9	4.5	95
IC2002-20	20G	Red	1.1	0.8	3.3	65
IC2002-22	22G	Blue	0.9	0.6	2.5	36
IC2002-24	24G	Yellow	0.7	0.5	1.9	23
IC2002-26	26G	Purple	0.6	0.45	1.9	17

Packaging: 50 pcs / box, sterile

Disposable Intravenous Catheter Without Port And Without Wings

(MDA. Reg. No: GB24437168617)



CATALOGUE NO.	SIZE	COLOR CODE	OD (mm)	ID (mm)	LENGTH (cm)	FLOW RATE (ml/min)
IC2003-14	14G	Orange	2.0	1.7	4.5	305
IC2003-16	16G	Grey	1.7	1.3	4.5	200
IC2003-18	18G	Green	1.3	0.9	4.5/3.2	90
IC2003-20	20G	Red	1.1	0.8	3.2	56
IC2003-22	22G	Blue	0.9	0.6	2.5	36
IC2003-24	24G	Yellow	0.7	0.5	1.9	23
IC2003-26	26G	Purple	0.6	0.45	1.9	17

Packaging: 50 pcs / box, sterile

Blood Giving Set with Airvent (Blood Transfusion / Intravenous Admin Set, Double Chamber)

(MDA. Reg. No: GB52945443217)



CATALOGUE NO.	TYPE
IC3050	Standard

Packaging: 100 sets / box, sterile

Blood Giving Set with Airvent (Blood Transfusion / Intravenous Admin Set)

(MDA. Reg. No: GB52945443217)



CATALOGUE NO.	COLOR
IC3052	Standard
Packaging: 100 sets / box, sterile	

Intravenous Admin Set with Airvent & Y-Connection with Protective Shield

(MDA. Reg. No: GB43549443017)



CATALOGUE NO.	SIZE
IC3053A-21	21G X 1 1/2"
IC3053A-23	23G X 1"
Packaging: 100 sets / box, sterile	

Intravenous Admin Set with Airvent & Latex Free Y-Connection

(MDA. Reg. No: GB43549443017)



CATALOGUE NO.	SIZE
IC3055-21	21G X 1 1/2"
IC3055-23	23G X 1"
Packaging: 100 Sets / box, sterile	

Intravenous Admin Set with Airvent Adult, 20 Drops Fitted with 21G Needle

(MDA. Reg. No: GB43549443017)

Optional:

- Available with luer-lock with rubber y-site



CATALOGUE NO.	SIZE
IC3056	21G X 1 1/2"
Packaging: 100 Sets / box, sterile	

Intravenous Admin Set with Airvent Paediatric, 60 Drops Fitted with 23G Needle

(MDA. Reg. No: GB43549443017)

Optional:

- Available with luer-lock with rubber y-site



CATALOGUE NO.	SIZE
IC3057	23G X 1"
Packaging: 100 Sets / box, sterile	

Microdrip Precision Solution Administration Set

(MDA. Reg. No: GB56356265017)



**DEHP
Free**

CATALOGUE NO.	SIZE
IC3058	110 ml
IC3059	150 ml
Packaging: 25 sets / box, sterile	

Peritoneal Dialysis Transfusion Set (Administration Set for Peritoneal Dialysis)

(MDA. Reg. No: GC490131028718)



CATALOGUE NO.	TYPE
IC3012-10	Standard
Packaging: 10 pcs / box, sterile	

Peritoneal Dialysis Catheter Set (Thocar Cannula for P.D.S.)

(MDA. Reg. No: GC490131028718)



CATALOGUE NO.	SIZE
IC3010-25A	Adult
IC3010-25C	Child
Packaging: 25 sets / box, sterile	

Extension Tube With Male Female Luer Lock

(MDA. Reg. No: GB77303168717)



CATALOGUE NO.	LENGTH
IC3036-200	200 cm
Packaging: 50 pcs / box, sterile	

CVP Manometer (MDA. Reg. No: GB953401028418)



CATALOGUE NO.	TYPE
IC3008	Standard
Packaging: Single pack, sterile	

Tube System For CVP Measurement Complete With Iv Admin Set
(MDA. Reg. No: GB953401028418)



CATALOGUE NO.	TYPE
IC3009	Standard
Packaging: Single pack, sterile	

Ryle's Tube (MDA. Reg. No: GB73663404917)

GASTROENTEROLOGY



CATALOGUE NO.	SIZE	COLOR CODE
IC4032-08	08 Fg	Blue
IC4032-10	10 Fg	Black
IC4032-12	12 Fg	White
IC4032-14	14 Fg	Green
IC4032-16	16 Fg	Orange
IC4032-18	18 Fg	Red
IC4032-20	20 Fg	Yellow
Packaging: 50 pcs / box, sterile		

Infant Feeding Tube (MDA. Reg. No: GB73663404917)



CATALOGUE NO.	SIZE	COLOR CODE
IC4013-05	5 Fg	Grey
IC4013-06	6 Fg	Green
IC4013-08	8 Fg	Blue
IC4013-10	10 Fg	Black
Packaging: 100 pcs / box, sterile		

Mucus Extractor (MDA. Reg. No: GB19674443117)



CATALOGUE NO.	SIZE
IC5018-08	8 Fg
IC5018-10	10 Fg
IC5018-12	12 Fg
IC5018-14	14 Fg
Packaging: Single pack, sterile	

Under Water Seal Pleural Drainage System (MDA. Reg. No: GB89133604018)



CATALOGUE NO.	SIZE
IC5020	2000 ml
IC5031	1200 ml
IC5032	500 ml
Packaging: Single pack, sterile	

Close Wound Suction Unit (MDA. Reg. No: GB18992603918)



CATALOGUE NO. (Capacity: 400-600 ml)	SIZE
IC5008-10	10 Fg
IC5008-12	12 Fg
IC5008-14	14 Fg
IC5008-16	16 Fg
IC5008-18	18 Fg
Packaging: Single pack, sterile	

Suction Connecting Tube with Adaptor (MDA. Reg. No: GB98020310517)



CATALOGUE NO.	LENGTH
IC551-020	2.0 m
IC551-030	3.0 m
Packaging: Single pack, sterile	


Yankaur Suction Set (MDA. Reg. No: GB789971130518)



CATALOGUE NO.	LENGTH
IC5016A	Standard
Packaging: Single pack, sterile	

Yankaur Suction Handle (Standard tip & Crown tip)

(MDA. Reg. No: GB789971130518)



CATALOGUE NO.	SIZE
IC9014DS	Standard tip
IC9014DC	Crown tip
Packaging: 25 pcs / box, sterile	

Colostomy Bag (MDA. Reg. No: GMD69269664018A)



CATALOGUE NO.	SIZE
IC4020	140mm x 250mm
Packaging: 30 pcs / pack, sterile	

Thoracic Trocar Catheter (MDA. Reg. No: GB97717265217)



CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC5030-08	8 Fg	IC5030-20	20 Fg
IC5030-10	10 Fg	IC5030-24	24 Fg
IC5030-12	12 Fg	IC5030-28	28 Fg
IC5030-14	14 Fg	IC5030-32	32 Fg
IC5030-16	16 Fg	IC5030-36	36 Fg
IC5030-18	18 Fg	IC5030-40	40 Fg
Packaging: Single pack, sterile			

Chest Drainage Catheter Straight (MDA. Reg. No: GB22318264417)



CATALOGUE NO.	SIZE
IC5010S-12	12 Fr
IC5010S-14	14 Fr
IC5010S-16	16 Fr
IC5010S-18	18 Fr
IC5010S-20	20 Fr
IC5010S-22	22 Fr
IC5010S-24	24 Fr
IC5010S-28	28 Fr
IC5010S-32	32 Fr
IC5010S-36	36 Fr
IC5010S-40	40 Fr
Packaging: Single pack, sterile	

Chest Drainage Catheter with Right Angle (MDA. Reg. No: GB22318264417)



CATALOGUE NO.	SIZE
IC5010A-12	12 Fr
IC5010A-16	16 Fr
IC5010A-20	20 Fr
IC5010A-24	24 Fr
IC5010A-28	28 Fr
IC5010A-32	32 Fr
Packaging: Single pack, sterile	

PATHOLOGY

Petri Dishes (MDA. Reg. No: GA41481135017)



Single compartment

Double compartment

CATALOGUE NO.	SIZE
IC24-3001R	90 mm x 15 mm (Single compartment)
IC24-3002	90 mm x 15 mm (Double compartment)
Packaging: 500 pcs / box	

Microscope Glass Slide

(MDA. Reg. No: IVD21875330317A)



CATALOGUE NO.	SIZE
IC24-2409R	75mm x 25mm (Ground Edge)
IC24-2409FBS	75mm x 25mm (One Frosted End at Both Side)
Packaging: 100 pcs / box	

Proctoscope

(MDA. Reg. No: GA8145919-27193)



* Patented

CATALOGUE NO.	SIZE	TYPE
IC8900A	Standard	Trunk End
Packaging: Single pack, sterile		

CATALOGUE NO.	SIZE	TYPE
IC8900	Standard	Fluted-beaked
Packaging: Single pack, non-sterile		

Polystyrene / Polypropylene Specimen Urine Container

(MDA. Reg. No: GMD24762217117A)



White

Yellow

Red

Specimen Urine Container (60ml)

CATALOGUE NO.	SIZE
IC1061RPP-NS (Urine Specimen Container)	60 ml
IC1061WPP-NS (Urine Specimen Container)	60 ml
IC1061YPP-NS (Urine Specimen Container)	60 ml
Packaging: 500 pcs / box, non-sterile	

CATALOGUE NO.	SIZE
IC1063RPS-S (Urine Specimen Container)	60 ml
IC1063WPS-S (Urine Specimen Container)	60 ml
IC1063YPS-S (Urine Specimen Container)	60 ml
Packaging: 500 pcs / box, sterile	

Polypropylene Plastic Container with Screw cap

(MDA. Reg. No: GMD4610983517A)



CATALOGUE NO.	SIZE
IC1062SC-90 (Urine Specimen Container)	90 ml
Packaging: 500 pcs, sterile	

Sputum Cup

(MDA. Reg. No: GMD93999690918A)



CATALOGUE NO.	SIZE
IC2054-500	2 oz
Packaging: 500 pcs / box	

Vaginal Speculum

(MDA. Reg. No: GMD69269664018A)



Small

Medium

Large

CATALOGUE NO.	SIZE
IC2043-S	Small
IC2043-M	Medium
IC2043-L	Large
Packaging: Single pack, sterile	

PROCEDURE & CLINICAL SETS

Basic Dressing Set (Premium)

(MDA. Reg. No: GMD3652322316A)



CATALOGUE NO.	TYPE
IC7000P108	3 Forceps
Packaging: Single pack, sterile	

Basic Dressing Set (Standard)

(MDA. Reg. No: GMD3652322316A)



CATALOGUE NO.	TYPE
IC7000S	2 Forceps
Packaging: Single pack, sterile	

Basic Dressing Set (Haemo)

(MDA. Reg. No: GMD3652322316A)



CATALOGUE NO.	TYPE
IC7000H	1 Forcep
Packaging: Single pack, sterile	

Vaginal Examination Set

(MDA. Reg. No: GMD1179279317A)



CATALOGUE NO.	TYPE
IC7001D	Loop
IC7001D-1	adhesive pad
Packaging: Single pack, sterile	

Catheterization Set

(MDA. Reg. No: GMD98384104017A)



CATALOGUE NO.	TYPE
IC7002D	Standard
Packaging: Single pack, sterile	

PROCEDURE & CLINICAL SETS

Delivery Set

(MDA. Reg. No: GMD15739209317A)



CATALOGUE NO.	TYPE
IC7003D	Standard
Packaging: Single pack, sterile	

Injection Set

(MDA. Reg. No: GMD4544579717A)



CATALOGUE NO.	TYPE
IC7005D	Standard
Packaging: Single pack, sterile	

Pap Smear Test Kit

(MDA. Reg. No: GMD244492217A)



CATALOGUE NO.	TYPE
IC7006D	Standard
Packaging: Single pack, sterile	

Toilet & Suture Set

(MDA. Reg. No: GMD25205167417A)



CATALOGUE NO.	TYPE
IC7007D	Standard
Packaging: Single pack, sterile	

Top & Tail / Baby Cleaning Set

(MDA. Reg. No: GMD9323693217A)



CATALOGUE NO.	TYPE
IC7008D	Standard
Packaging: Single pack, sterile	

Eye Dressing Set

(MDA. Reg. No: GMD48257774918A)



CATALOGUE NO.	TYPE
IC7010D	Standard
Packaging: Single pack, sterile	

Cord Dressing Set (MDA. Reg. No: GMD22439210017A)



CATALOGUE NO.	TYPE
IC7011D-3	Standard
Packaging: Single pack, sterile	

Suction Set (MDA. Reg. No: GMD44504260217A)



CATALOGUE NO.	SIZE
IC7013D	Standard
Packaging: Single pack, sterile	

Blood Taking Set (MDA. Reg. No: GMD31881210117A)



CATALOGUE NO.	TYPE
IC7020	Standard
Packaging: Single pack, sterile	

Swabbing Set (MDA. Reg. No: GMD72303105317A)



CATALOGUE NO.	TYPE
IC7021B	Standard
Packaging: Single pack, sterile	

PROCEDURE SET COMPONENTS

Kidney Dish (MDA. Reg. No: GMD77454785318A)



CATALOGUE NO.	SIZE	COLOUR
IC2020-PB	700 ml	BLUE
IC2020-PG	700 ml	GREEN
Packaging: Single pack		

Forceps (MDA. Reg. No: GMD28172787718A)



CATALOGUE NO.	SIZE	COLOUR
IC6029B	12.5 cm	BLUE
IC8029R	12.5 cm	RED
Packaging: 100 pcs / pack		

Umbilical Cord Clamp (MDA. Reg. No: GA5125083616)



CATALOGUE NO.	SIZE
IC5084	Standard
Packaging: 100 pcs / box, sterile	

WOUND CARE & DRESSINGS

Cotton Wool Roll (MDA. Reg. No: GMD16545789218A)



CATALOGUE NO.	SIZE
IC2158-500	500 grams
Packaging: single pack	

Combine Dressing Roll-Gamgee (MDA. Reg. No: GMD76231428818A)



CATALOGUE NO.	SIZE
IC2159-500	500 grams (20 cm x 10 m)
Packaging: single pack	

Gauze Ribbon (MDA. Reg. No: GMD77112487418A)



CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC2093-122024	1.25 cm x 5.5 m	IC2093-902024	5 cm x 90 m
IC2093-551915	5 cm x 5.5 m	IC2093-102024	2.5 cm x 100 m
Packaging: single pack			

Lint Absorbent (MDA. Reg. No: GMD44226404918A)



CATALOGUE NO.	SIZE
IC2180-500	500 grams
Packaging: single pack	

Absorbent Gauze Swab (MDA. Reg. No: GMD77112487418A)



CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC2058-501915	5 cm x 5 cm-8 Ply	IC2068-751915	7.5 cm x 7.5 cm-32 Ply
IC2058-751915	7.5 cm x 7.5 cm-8 Ply	IC2063-752024	7.5 cm x 7.5 cm-24 Ply
IC2063-751915	7.5 cm x 7.5 cm-24 Ply	IC2067-752024	7.5 cm x 7.5 cm-32 Ply
Packaging: 100 pcs / pack			

Plaster Extension (MDA. Reg. No: GMD31219498518A)



CATALOGUE NO.	SIZE
IC093003-02	10 cm x 2.5 m
Packaging: single pack	

Skin Color Crepe Bandage (MDA. Reg. No: GMD24431526418A)



CATALOGUE NO.	SIZE
IC2121-050	5 cm x 4 M
IC2121-075	7.5 cm x 4 M
IC2121-100	10.0 cm x 4 M
IC2121-150	15.0 cm x 4 M

Packaging: 12 rolls / pack

Dental Wool Roll No. 2 (MDA. Reg. No: GMD76257427118A)



CATALOGUE NO.	SIZE
IC2152	10 mm x 38 mm

Packaging: 500 pcs / box

Cotton Wool Balls (MDA. Reg. No: GMD16545789218A)



CATALOGUE NO.	SIZE
IC2151-005A	0.5 grams

Packaging: 250 pcs / pack

Transparent Film Dressing (MDA. Reg. No: GMD45610109117A)



CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC2162-010	10 cm x 12 cm	IC2162-065	6.5 cm x 5 cm
IC2162-010B	10 cm x 20 cm	IC2162-095	9.5 cm x 8.5 cm
IC2162-010C	10 cm x 25 cm	IC2162-155	15.5 cm x 8.5 cm
IC2162-040	4 cm x 5 cm	IC2162-235	23.5 cm x 8.5 cm
IC2162-060	6 cm x 7 cm	Packaging: 100 pcs / box, sterile	

Transparent Infusion Set Dressing (MDA. Reg. No: GMD42126130817A)



CATALOGUE NO.	SIZE
IC2163-067	6 cm x 7 cm

Packaging: 100 pcs / box, sterile

Eye Pad (MDA. Reg. No: GMD1927680617A)



Standard

Adhesive

CATALOGUE NO. (STANDARD)	CATALOGUE NO. (ADHESIVE)	SIZE
IC2150A	IC2149A	7 cm x 9 cm (Adult)
IC2150C	IC2149C	5 cm x 6.2 cm (Child)

Packaging: 50 pcs / box, sterile

WOUND CARE & DRESSINGS

Paraffin Gauze Dressing (MDA. Reg. No: GMD62261109417A)



CATALOGUE NO.	SIZE
IC5.02	10 cm x 10 cm
Packaging: 10 pcs / box, sterile	

Soft Paper Surgical Tape (MDA. Reg. No: GMD45066451418A)



CATALOGUE NO. (Soft Paper Surgical Tape)	SIZE	PACKAGING
IC2110-125	1.25 cm x 9.1 m	24 pcs / box
IC2110-250	2.5 cm x 9.1 m	12 pcs / box
IC2110-500	5.0 cm x 9.1 m	6 pcs / box
IC2110-750	7.5 cm x 9.1 m	4 pcs / box

4 Folded BPC Gauze Rolls (MDA. Reg. No: GMD76231428818A)



CATALOGUE NO.	SIZE
IC2163-1915	90 cm x 100 Y
IC2163-2024	90 cm x 100 Y
IC2163-XRD	90 cm x 100 Y
Packaging: single pack, sterile	

Orthopaedic Bandage (MDA. Reg. No: GMD24431526418A)



CATALOGUE NO.	SIZE
IC2157-050	5 cm x 2.75 M
IC2157-075	7.5 cm x 2.75 M
IC2157-100	10 cm x 2.75 M
IC2157-150	15 cm x 2.75 M
IC2157-200	20 cm x 2.75 M
Packaging: 10 pcs / pack	

Male Urinal Container, Plastic (with cover) (MDA. Reg. No: GMD12875361617A)



CATALOGUE NO.	SIZE
IC2046MC	1000 ML
Packaging: single pack	

Bed Pan, Plastic (MDA. Reg. No: GMD12875361617A)



CATALOGUE NO.	SIZE
IC2045	Standard
Packaging: single pack	

Alcohol Swabs

(MDA. Reg. No: GMD27041357317A)



CATALOGUE NO.	SIZE
IC0645	2 ply-30 x 30 mm (folded)
Packaging: 200 pcs / box	

PHARMACEUTICAL PLASTIC CONTAINER

Medicine Cup

(MDA. Reg. No: GMD39392630718A)



CATALOGUE NO.	SIZE
MC002-20ML	20 ml
Packaging: 100 pcs / pack	

Lotion Bottle



CATALOGUE NO.	SIZE	PACKAGING
IC9005-060	60 ml	500 pcs / pack
IC9005-090	90 ml	
IC9005-120	120 ml	250 pcs / pack
IC9005-240	240 ml	

Poison Bottle



CATALOGUE NO.	SIZE	PACKAGING
IC9003-060	60 ml	500 pcs / pack
IC9003-090	90 ml	
IC9003-120	120 ml	250 pcs / pack
IC9003-240	240 ml	

Ointment Box



CATALOGUE NO.	SIZE	PACKAGING
IC9004-015	1/2 oz (15 g)	144 pcs / pack
IC9004-030	1 oz (30 g)	144 pcs / pack
IC9004-120	4 oz (120 g)	144 pcs / pack
IC9004-240	8 oz (240 g)	Single Pack
IC9004-480	16 oz (480 g)	Single Pack

Mustajab Bottle with Flip Top Cap



CATALOGUE NO.	SIZE	PACKAGING
IC9007-100	100 ml	500 pcs / pack
IC9007-130	130 ml	
IC9007-200	200 ml	250 pcs / pack
IC9007-260	260 ml	

Eye / Ear Drop (MDA. Reg. No: GMD16117785218A)



CATALOGUE NO.	SIZE
IC9009-005	5 ml
IC9009-010	10ml
Packaging: 144 pcs / pack	

Dispensing Bottle



CYLINDRICAL			RECTANGLE		
CATALOGUE NO.	SIZE	PACKAGING	CATALOGUE NO.	SIZE	PACKAGING
IC9002C-030	30 ml	1000 pcs / pack	IC9002R-030	30 ml	1000 pcs / pack
IC9002C-060	60 ml	500 pcs / pack	IC9002R-060	60 ml	500 pcs / pack
-	-	-	IC9002R-090	90 ml	500 pcs / pack
IC9002C-120	120 ml	500 pcs / pack	IC9002R-120	120 ml	500 pcs / pack
IC9002C-240	240 ml	250 pcs / pack			

Also available in:

	Red cap
	Blue cap
	White cap

Wheelchair



Available in:
Solid & Pneumatic Tyre

REHABILITATION

CAT. NO.	SIZE
IC869X-46	Standard
Packaging: Per unit	


Underpad



CATALOGUE NO.	SIZE	PACKAGING
IC2047B-150	30"x30"	150 pcs / carton
IC2047C-300	17"x24"	300 pcs / carton

OTHERS

Surgical Disposable Face Mask (3-ply) (MDA. Reg. No: GMD96851528118A)



Tie-on-type



Earloop-type

CATALOGUE NO.	TYPE	COLOR
IC2102-003TG	3-ply (tie-on-type)	green
IC2102-003EG	3-ply (earloop-type)	green
Packaging: 50 pcs / box		

Surgical Disposable Face Mask (non-woven fabric) (MDA. Reg. No: GMD96851528118A)



Tie-on-type

CATALOGUE NO.	TYPE	COLOR
IC2102-003TGW	3-ply (tie-on-type)	green
Packaging: 50 pcs / box		

Surgical Disposable Face Mask (kids) (MDA. Reg. No: GMD96851528118A)



Earloop-type

CATALOGUE NO.	TYPE	COLOR
IC2102-003EGP	3-ply (earloop-type/kids)	green
Packaging: 50 pcs / box		

Surgical Disposable Face Mask (4-ply) (MDA. Reg. No: GMD96851528118A)



Tie-on-type



Earloop-type

CATALOGUE NO.	TYPE	COLOR
IC2017-004TB	4-ply (tie-on-type)	black
IC2017-004EL	4-ply (earloop-type)	black
Packaging: 50 pcs / box		

IDEALCARE N95/FFP2 MASK (Disposable Particulate Respirator)

(MDA. Reg. No: GA3713521-65410)



Earloop



Headloop (adjustable)

CATALOGUE NO.	SIZE
IC2107-EL	Earloop
IC2107-HL	Headloop (adjustable)
Packaging: 50 pcs / box	

Medical Face Shield / Isolation Mask (MDA. Reg. No: GA5277621-68822)


CATALOGUE NO.	SIZE
IC0032	22cm x 33cm
Packaging: single pack	

PPE Protective Coverall Suit (MDA. Reg. No: GA5937921-68882)
Protective
coverall suit

Boot cover

CATALOGUE NO.	SIZE
IC2104-L (Protective coverall suit)	170cm
IC2104-XXL (Protective coverall suit)	180cm
Packaging: single pack, sterile	

CATALOGUE NO.	SIZE
IC2106	standard
Packaging: single pack, sterile	

Disposable Protective Coverall (MDA. Reg. No: GA5937921-68882)



CATALOGUE NO.	SIZE
IC2108-165	M
IC2108-170	L
IC2108-175	XL
IC2108-180	XXL
IC2108-185	XXXL

Packaging: single pack

Disposable Boot Cover (MDA. Reg. No: GA5670421-68824)



CATALOGUE NO.	TYPE	SIZE
IC2100-000B	42gsm	standard

Packaging: 50 pcs (25 pairs) / pack

Disposable Shoe Cover (MDA. Reg. No: GA4661021-66726)



CATALOGUE NO.	TYPE	SIZE
IC2105-000B	42gsm	standard

Packaging: 100 pcs / pack

Disposable Head Cover (Surgical Cap / Nursing Cap)



CATALOGUE NO.	TYPE	SIZE
IC2095-030B	30gsm	standard
IC2095-000B	42gsm	standard

Packaging: 100 pcs / pack

Disposable Hood Cover (Tudung Cover) (MDA. Reg. No: GA4326721-56820)


CATALOGUE NO.	SIZE
IC2094-000B	standard
Packaging: 50 pcs / pack	

Disposable Plastic Apron (MDA. Reg. No: GA4741521-67936)


CATALOGUE NO.	SIZE
IC2008-50M	With sleeve and thumb hook
IC2008-140	Without sleeve
Packaging: single pack	

Non-Woven Isolation Gown (MDA. Reg. No: GA4260519-27772)


CATALOGUE NO.	MATERIAL	SIZE
IC2097-NWBS	SMS	125x145cm / blue
IC2097-NWG	SMMS	125x145cm / blue
IC2097-NWB	PP	115x130cm / blue
Packaging: single pack		

AUTHORIZE REPRESENTATIVE
(POLY MEDICURE LIMITED)

Guedel Airway Oro-Pharyngeal Airway (MDA. Reg. No: GB13993120217)



CATALOGUE NO.	SIZE	COLOR CODE
20050	000 (40mm)	Purple
20051	00 (50mm)	Blue
20052	0 (60mm)	Black
20053	1 (70mm)	White
20054	2 (80mm)	Green
20055	3 (90mm)	Orange
20056	4 (100mm)	Red
20057	5 (110mm)	Yellow

Packaging: 50 pcs / box, sterile

Catheter Mount (MDA. Reg. No: GB71880264617)



CATALOGUE NO.	SIZE
20180	I.D 15 mm / O.D 22 mm

Packaging: single pack, sterile

Spinal Needle (MDA. Reg. No: GD21209204317)



CATALOGUE NO. (Spinal Needle with Whitacre Point Needle)	CATALOGUE NO. (Spinal Needle with Quincke Point Needle)	SIZE	CATALOGUE NO. (Spinal Needle with Whitacre Point Needle)	CATALOGUE NO. (Spinal Needle with Quincke Point Needle)	SIZE
20400	20919	16 G	20406	20925	22 G
20401	20920	17 G	20407	20926	23 G
20402	20921	18 G	20408	20927	24 G
20403	20922	19 G	20409	20928	25 G
20404	20923	20 G	20917	20929	26 G
20405	20924	21 G	20918	20930	27 G

Packaging: Single pack, sterile

INFUSION & TRANSFUSION

Polyway Pro Click-Three Way Stop Lipid Resistant With Stop Cock
Resistant With 45 Degree Click (MDA. Reg. No: GB15021604118)



CATALOGUE NO.	COLOR
13058B	blue
13058R	red

Packaging: 50 pcs / box, sterile

Polyspike (MDA. Reg. No: GA10447219-34401)



Polyspike Types & Codes		Fluid Filter	Air Filter
	13490	5µ	0.1µ
	13491	5µ	0.1µ
	13068	NA	0.1µ

CATALOGUE NO.	SIZE
13490	Standard Cap
13491	Standard Cap
13068	Standard Cap
Packaging: 100 pcs / box	

Extension Tubes / Lines

MDA. Reg. No: GB69314264517)



VARIANTS				SIZE
CATALOGUE NO. (PVC High Pressure Extension Line, ID - 1.0 mm x OD - 3.0 mm)	CATALOGUE NO. (PVC Low Pressure Extension Line, ID - 3.0 mm x OD - 4.1 mm)	CATALOGUE NO. (PVC Free (PE) High Pressure Extension Line (straight), ID - 1.0 mm x OD - 2.0 mm)	CATALOGUE NO. (PVC Free (PE) High Pressure Extension Line (coiled), ID - 1.0 mm x OD - 2.0 mm)	
13120	-	13125	-	25 cm
13121	13183	13126	-	50 cm
13122	13185	13127	13130	100 cm
13123	13187	13128	-	150 cm
13124	13189	13129	13131	200 cm
-	-	-	13132	300 cm
Packaging: 50 pcs / box, sterile				

Polyway with Extention Tube-Three Way Stop Cock With Extension Tube

(MDA. Reg. No: GB69314264517)



CATALOGUE NO.	SIZE
13016	10 cm
13017	25 cm
Packaging: 75 pcs / box, sterile	
13018	50 cm
Packaging: 50 pcs / box, sterile	
13020	100 cm
13021	150 cm
13022	200 cm
Packaging: 40 pcs / box, sterile	

Oncofusion Set (IV Infusion Set with 0.2 Micron In-Line and Priming Filter)

(MDA. Reg. No: GB88419265117)



CATALOGUE NO.	SIZE
14163	150 cm
Packaging: Single pack, sterile	

Autofusion Set (I.V. Infusion Set with Auto Air Stop and Priming Filter)

(MDA. Reg. No: GB88419265117)



CATALOGUE NO.	SIZE
14261	180 cm
Packaging: Single pack, sterile	

Polytrol-Micro Set (I.V. Infusion Set with Micro Drip and Flow Regulator)

(MDA. Reg. No: GB88419265117)



CATALOGUE NO.	SIZE
14261	150 cm
Packaging: Single pack, sterile	

Photofusion Set (I.V. Infusion Set for Photosensitive Drugs)

(MDA. Reg. No: GB88419265117)



CATALOGUE NO.	SIZE
14230	150 cm
Packaging: Single pack, sterile	

Safe Conekt (MDA. Reg. No: GB2186920-43643)



CATALOGUE NO.	SIZE
13326	Large bore Mono-10cm
13323	Large bore Duo-15cm
13321	Small bore Duo-13cm
13322	Small bore Trio-15cm
13324	Large bore Trio-15cm
Packaging: Single pack, sterile	

Injection Stopper (MDA. Reg. No: GMD3146469216A)



CATALOGUE NO.	SIZE
13354	Standard
Packaging: 300 pcs / box, sterile	

Insulin Syringe (Blister Pack)-U/100 (MDA. Reg. No: GB834271685517)



CATALOGUE NO.	SIZE
51040	29 G
51050	30 G
51060	31 G
Packaging: 50 pcs/box, sterile	

A.V. Fistula Needle (MDA. Reg. No: GB71204168417)



DIALYSIS

SIZE	CATALOGUE NO. (Fixed Wings)	CATALOGUE NO. (Rotating Wings)
15 G	90030	90230
16 G	90031	90231
17 G	90032	90232
Packaging: Single pack, sterile		

Levin's Tube-Gastro-Duodenal Feeding Tube

(MDA. Reg. No: GB52652168317)




GASTROENTEROLOGY

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
40301	8 Fg	40306	18 Fg
40302	8 Fg	40307	18 Fg
40303	10 Fg	40308	20 Fg
40304	12 Fg	40309	22 Fg
40305	14 Fg	40310	24 Fg
Packaging: 50 pcs / box, sterile			

High Pressure Vacuum Bottle

(MDA. Reg. No: GMD8443663216A)



SURGERY & WOUND DRAINAGE

CATALOGUE NO.	SIZE
90510	200 ml
90511	400 ml
90512	600 ml
Packaging: Single pack, sterile	

CENTRAL VENOUS ACCESS CATHETER

Central Venous Catheter (MDA. Reg. No: GD28415264917)

CATALOGUE NO.	SIZE	LENGTH	VARIANTS	
15001	3 Fr / 20 G	80 mm	SINGLE LUMEN	Packaging: Single pack, sterile
15002	3 Fr / 20 G	130 mm		
15003	3 Fr / 20 G	150 mm		
15004	3 Fr / 20 G	160 mm		
15005	3 Fr / 20 G	200 mm		
15006	4 Fr / 18 G	130 mm		
15007	4 Fr / 18 G	150 mm		
15008	4 Fr / 18 G	160 mm		
15009	4 Fr / 18 G	200 mm		
15010	5 Fr / 16 G	160 mm		
15011	5 Fr / 16 G	200 mm		
15012	5 Fr / 16 G	300 mm		
15013	6 Fr / 14 G	160 mm		
15014	6 Fr / 14 G	200 mm		
15015	6 Fr / 14 G	300 mm		
15016	7 Fr / 12 G	160 mm		
15017	7 Fr / 12 G	200 mm		
15018	7 Fr / 12 G	300 mm		
15019	8 Fr / 12 G	200 mm		
15020	9 Fr / 12 G	200 mm		
15050	4 Fr / 22 G(P) / 22 G(D)	130 mm	DOUBLE LUMEN	
15051	4 Fr / 22 G(P) / 22 G(D)	160 mm		
15052	4 Fr / 22 G(P) / 22 G(D)	200 mm		
15053	5 Fr / 18 G(P) / 20 G(D)	130 mm		
15054	5 Fr / 18 G(P) / 20 G(D)	200 mm		
15055	7 Fr / 16 G(P) / 16 G(D)	160 mm		
15056	7 Fr / 16 G(P) / 16 G(D)	200 mm		
15057	7 Fr / 16 G(P) / 16 G(D)	300 mm		
15058	7 Fr / 14 G(P) / 18 G(D)	160 mm		
15059	7 Fr / 14 G(P) / 18 G(D)	200 mm		
15060	7 Fr / 14 G(P) / 18 G(D)	300 mm		
15061	8 Fr / 14 G(P) / 14 G(D)	160 mm		
15062	8 Fr / 14 G(P) / 14 G(D)	200 mm		
15063	8 Fr / 14 G(P) / 14 G(D)	300 mm		
15100	5.5 Fr / 22 G(P) / 20 G(D) / 22 G(M)	80 mm	TRIPLE LUMEN	
15101	5.5 Fr / 22 G(P) / 20 G(D) / 22 G(M)	160 mm		
15102	5.5 Fr / 22 G(P) / 20 G(D) / 22 G(M)	200 mm		
15103	7 Fr / 18 G(P) / 16 G(D) / 18 G(M)	160 mm		
15104	7 Fr / 18 G(P) / 16 G(D) / 18 G(M)	200 mm		
15105	7 Fr / 18 G(P) / 16 G(D) / 18 G(M)	300 mm		
15106	4.5 Fr / 22 G(P) / 20 G(D) / 22 G(M)	60 mm		
15107	12 Fr / 12 G(P) / 16 G(D) / 12 G(M)	150 mm		
15564	7 Fr / 18 G(P) / 16 G(D) / 18 G(M)	200 mm	TRIPLE LUMEN NO VACENT TRIO PRO-KINK RESISTANT GUIDE WIRE	



VISION

To be a one-stop centre for hospitals, healthcare centers, clinics and pharmacies. Already a market leader in Malaysia and with the continuous improvement attitude it will become a market leader regionally as well as internationally.

MISSION

Support and enhance the market leadership of the Disposable Product Businesses through the distribution of the highest quality products at superior service levels and most competitive prices.

The foundation for this effort will be the full partnership of our employee in the attainment of our goals through open communications and total involvement in our business.



Contact Us



Northern Territory

Business Development Department

Tel.: 011-6503 5775

Email: sales01@idealcare.com.my

East Coast Territory

Business Development Department

Tel.: 011-6308 5775

Email: sales04@idealcare.com.my

Central & Southern Territory

Business Development Department

Tel.: 011-6343 5775

Email: marketing@idealcare.com.my

East Malaysia (Sabah/Sarawak)

Business Development Department

Tel.: 011-6321 5775

Email: sales02@idealcare.com.my

Customer Care

Business Development Department

Tel.: 019-415 7775

Email: customercare@idealcare.com.my

Ideal Health Care Sdn. Bhd. 1999901009961 (484861-V)

Plant 1



• Distribution Centre:

No. 70 & 71, Jalan Sungai Tukang 2/1,
Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 422 7775, (6) 04 - 422 8775
Fax : 6 (04) - 422 8775
E-mail : info@idealcare.com.my

Plant 2



• Head Office / Center of Excellence, COE

Lot PT. 80598,
Kawasan Perindustrian LPK Fasa 1,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 444 6875, (6) 04 - 444 6975
Fax : 6 (04) - 444 5975
E-mail : info@idealcare.com.my



www.idealcare.com.my