

INFRARED FOREHEAD THERMOMETER

PPE COVID 19

- ✓ One second quick test.
- ✓ The error is less than $\pm 0.2^{\circ}\text{C}$.
- ✓ Professional infrared core.
- ✓ Packaging: single pack.

Technical Specs	
Measure range:	32°C - 42.9°C
Battery:	DC 3.0V
Electricity prompt:	2.4V \pm 0.2V
Memory:	32 sets
Backlight:	3 colors
Measure time:	one second
Auto power-off:	8s



IHSB REF. NO.

IC0646-GW100

SIZE

Standard

Ideal Healthcare Sdn. Bhd. (494861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 422 7775 / 422 8775 Fax : 6 (04) - 422 6775
Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

"Your One Stop Centre for Medical & Healthcare Products"

www.idealcare.com.my

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 13485 =



IHSB-CAT-BDD-PPE-06-00

