

**IDEALCARE N95/FFP2 MASK
(DISPOSABLE PARTICULATE RESPIRATOR)**

SURGICAL RESPIRATOR

Common Features:

- ▶ Used as respiratory protective device designed to achieve a very close facial fit and very efficient filtration of airborne particles meant for medical or non-medical use
- ▶ Disposable and single use reduce the risk recontamination and cross-infection
- ▶ The product is fluid resistant respirator and tested by National Institute of Occupational Safety and Health (NIOSH) and cleared by MDA Malaysia as medical device
- a) Air-purifying medical respirators with filter efficiency level above 98%**
- b) EN149 (LEVEL FFP2)/ASTM F2100-19 compliance**
- c) Medical grade N95/FFP2 mask**
- ▶ Non-valve respirator
- ▶ Cushioned foam/sponge around the nose clip gives maximum user comfort
- ▶ An adjustable nose clip provides custom fit
- ▶ Lightweight construction facilitates a long wear time
- ▶ Advanced electrostatic media is designed for ease of breathing
- ▶ Style: Duck-Bill
- ▶ Material is from soft lining and non-irritating to the skin
- ▶ **STERILE** and individual pack provides maximum safety and prevent contamination during storage
- ▶ To be worn by healthcare professionals during surgical procedure



ANTI DUST



ANTI DROPLETS



COMFORTABLE BREATHING



3 DIMENSIONAL STRUCTURE



SEAMLESS CRIMPING





PREMIUM MATERIAL





5 LAYER



ADJUSTABLE NOSE BRIDGE



• Earloop type



• Headloop type (Adjustable)
(Two-strap design with ultrasonic welded dual point attachment provides a secure seal)



Cushioned foam/sponge around the nose clip gives maximum user comfort



5 LAYER DESIGN



MDA Reg. No: GMD96851528118A

STERILE

IHSB REF. NO

IC2107-EL

IC2107-HL

TYPE

Earloop

Headloop (Adjustable)

Packaging: Single Pack, Sterile

SKU: 50 Pcs / Box



Ideal Healthcare Sdn. Bhd.(484861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) - 422 7775 / 422 8775 Fax : 6 (04) - 422 6775
Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

"Your One Stop Centre for Medical & Healthcare Products"

www.idealcare.com.my



idealcareMY

IHSB-CAT-BDD-PLY-02-00