

RESPIRATORY CARE

Oxygen Mask with Tubing (IC2020A, IC2020P, IC2020I & IC2020N)

- Designed for convenient oxygen therapy.
- Moulded face mask is made of medical grade PVC.
- Elongated mask for long term use.
- Clear and soft mask for patient comfort.
- Adjustable elastic strip and intergrated nose clip for proper positioning of mask.



Oxygen Mask High Concentration Non Rebreathing Mask (IC110A & IC1112P)

- Equip with soft, transparent and odour free mask.
- Light weight reservoir improves the oxygen efficiency.
- Suitable for high concentration oxygen therapy.
- Comes with a transparent plastic reservoir attached to the mask allows deep breathing without wastage of oxygen.




Nebulizer Mask with Chamber & Tubing (IC1120, IC1122 & IC1123)

- Soft and clear with anatomical form is ideal for long term use.
- Intergrated nose bridge for extra comfort.
- Adjustable elastic strip and intergrated nose clip for proper positioning of mask.
- Complete kit consists of : nebulizer mask, nebulizer chamber and tubing.



Nasal Oxygen Cannula (IC2016A, 2016P, 2016I & 2016N) (IC2017A, 2017P, 2017I & 2017N)

- Easy application and efficient administration of oxygen for patient's comfort.
- Suitable for easy and efficient administration of oxygen.
- Soft and kink resistant PVC tubing.
- Soft twin prong nasal tips to ensure equal distribution of oxygen through both air passages which provides maximum patient comfort.



* for illustration purpose only



Ideal Healthcare Sdn. Bhd.(484861-V)
 No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
 08000 Sungai Petani, Kedah, Malaysia.
 Tel : 6 (04) - 422 7775 / 422 8775 Fax : 6 (04) - 422 6775
 Email : info@idealcare.com.my

E. L. NO.: MDA-0003-KD314

www.idealcare.com.my

COMPANY WITH
 QUALITY SYSTEM
 CERTIFIED BY DNV GL
 = ISO 13485 =








idealcareMY

IHSB-CAT-BDD-RES-01-00

Oxygen Air Entrainment (IC7411& 7412)

- Made of medical grade PVC.
- Adjustable nose clip.
- Easy adjustable, variable oxygen air entrainment system ensure fixed and accurate concentrations of oxygen.
- Color coded, air entrainment low and medium concentration diluters.

Adjusted Oxygen Concentration	Color Code
24% (Blue)	
28% (Yellow)	
35% (Dark Green)	
40% (Pink)	
60% (Red)	



Oxygen Mask with Tubing

CATALOGUE NO.	SIZE
IC2020A	Elongated Adult
IC2020P	Elongated Paediatric
IC2020I	Elongated Infant
IC2020N	Elongated Neonate

Packaging: Single pack, 100 pcs / carton

Oxygen Mask High Concentration Non Rebreathing Mask

CATALOGUE NO.	SIZE
IC1110A	Adult
IC1112P	Paediatric

Packaging: Single pack, 100 pcs / carton

Nebulizer Mask with Chamber & Tubing

CATALOGUE NO.	SIZE
IC1120	Elongated Adult
IC1122	Elongated Paediatric
IC1123	Standard Neonate

Packaging: Single pack, 100 pcs / carton

Oxygen Air Entrainment

CATALOGUE NO.	SIZE
IC7411	Adult
IC7412	Paediatric

Packaging: Single pack, 100 pcs / carton

Nasal Oxygen Cannula

CATALOGUE NO.	SIZE	CATALOGUE NO.	SIZE
IC2016A	Adult	IC2017A	Adult
IC2016P	Paediatric	IC2017P	Paediatric
IC2016I	Infant	IC2017I	Infant
IC2016N	Neonate	IC2017N	Neonate

Packaging: Single pack, 100 pcs / carton **Packaging: Single pack, sterile, 100 pcs / carton**