

4 PLY SURGICAL FACE MASK

Components

1. Nose bar.
2. Tie-on-type & Earloop-type.
3. Ultrasonic seam.
4. Layers of medical gauze & bacterial filter.
5. Imprinted **IDEALCARE®** logo.



General Specifications

- Meets EN 14683 Type IIR and / ASTM F2100-11 LEVEL 3 standard.
- Bacterial Filtration Efficiency (BFE) ≥98%.
- Manufactured in ISO 7 clean room environment.
- Conforms and passes Bioburden testing.
- Maximum barrier protection for any situation with potential of exposure to heavy levels of aerosols, sprays and fluids.
- Provided excellent breathability.
- Assembled with exceptional filtration layer.
- Secured with maximum comfort.

New

***Also available in Earloop.**

Manufacturer Information

- IDEAL HEALTHCARE SDN. BHD.**
- Establishment License No: MDA-1080-KD120.
 - MDA Registration No.: **GMD96851528118A.**
 - Shelf Life: 5 years from the manufacturing date.

Quality Certification



Black

Product Specification

Type	TIE-ON	EARLOOP
Product code	IC2017-004TB	IC2017-004EL
Illustration for application		
Packaging	50 pcs / box X 40 boxes / carton	
Packaging illustration		



Ideal Healthcare Sdn. Bhd. (484861-V)
No. 70 & 71, Jalan Sungai Tukang 2/1, Kawasan Perusahaan Sungai Tukang,
08000 Sungai Petani, Kedah, Malaysia.
Tel : 6 (04) – 422 7775 / 422 8775 Fax : 6 (04) – 422 6775
Email : info@idealcare.com.my

E. L. NO.: MDA-1080-KD120

"Your One Stop Centre for Medical & Healthcare Products"

www.idealcare.com.my

